

Product Datasheet

Human FGF9 ELISA Kit (orb406124)

Description	Human fibroblast growth factor 9 (glia-activating factor) (FGF9) ELISA kit
Reactivity	Human
Range	6.25 pg/mL-400 pg/mL
Target	FGF9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P31371
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.