

## Product Datasheet

### Human estrogen ELISA Kit (orb406132)

|                          |                              |
|--------------------------|------------------------------|
| <b>Description</b>       | Human estrogen,E ELISA Kit   |
| <b>Reactivity</b>        | Human                        |
| <b>Range</b>             | 25 pg/mL-1000 pg/mL          |
| <b>Target</b>            | E                            |
| <b>Note</b>              | For research use only        |
| <b>Application notes</b> | Detection Wavelength: 450 nm |
| <b>Sample Types</b>      | serum, plasma                |
| <b>Sensitivity</b>       | 25 pg/mL                     |
| <b>Expiration Date</b>   | Please enquire.              |