

## Product Datasheet

### Human CCSA-2 ELISA Kit (orb406150)

|                          |                                                            |
|--------------------------|------------------------------------------------------------|
| <b>Description</b>       | Human colon cancer-specific antigen-2,CCSA-2 ELISA kit     |
| <b>Reactivity</b>        | Human                                                      |
| <b>Range</b>             | 3.33 ng/mL-80 ng/mL                                        |
| <b>Target</b>            | CCSA-2                                                     |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates |
| <b>Sensitivity</b>       | 1.33 ng/mL                                                 |
| <b>Expiration Date</b>   | Please enquire.                                            |