

Product Datasheet

Human CYP7A1 ELISA Kit (orb406159)

Description	Human Cholesterol 7 α -Hydroxylase (CYP7A)ELISA Kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Target	CYP7A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P22680
Sensitivity	4.7 pg/mL
Expiration Date	Please enquire.