

Product Datasheet

Human FAM132A ELISA Kit (orb406170)

Description	Human Protein FAM132A(FAM132A) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	FAM132A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q5T7M4
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.