

Product Datasheet

Human Dermcidin ELISA Kit (orb406184)

Description	Human Proteolysis Inducing Factor/Dermcidin(PIF/DCD)ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	DCD
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, tissue homogenates
Uniprot ID	P81605
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.