

## Product Datasheet

### Human Versican ELISA Kit (orb406218)

|                          |                                         |
|--------------------------|-----------------------------------------|
| <b>Description</b>       | human versican/PG-M/PG-350<br>ELISA kit |
| <b>Reactivity</b>        | Human                                   |
| <b>Range</b>             | 0.04 ng/ml-30 ng/ml.                    |
| <b>Target</b>            | VCAN                                    |
| <b>Note</b>              | For research use only                   |
| <b>Application notes</b> | Detection Wavelength: 450 nm            |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates.   |
| <b>Uniprot ID</b>        | <b>P13611</b>                           |
| <b>Sensitivity</b>       | 0.04 ng/ml.                             |
| <b>Expiration Date</b>   | Please enquire.                         |