

Product Datasheet

Human SFTPA1 ELISA Kit (orb406247)

Description	Human Pulmonary surfactant-associated protein A,SP-A ELISA Kit
Reactivity	Human
Range	15.6 pg/mL-1000pg/mL
Target	SFTPA1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine
Uniprot ID	Q8IWL2
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.