

Product Datasheet

Human SCGB3A2 ELISA Kit (orb406257)

Description	Human Secretoglobin family 3A member 2(SCGB3A2) ELISA kit
Reactivity	Human
Range	46.88 pg/mL-3000 pg/mL
Target	SCGB3A2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q96PL1
Sensitivity	11.72 pg/mL
Expiration Date	Please enquire.