

Product Datasheet

Human SFRP1 ELISA Kit (orb406258)

Description	Human Secreted Frizzled-Related Protein 1(SFRP1) ELISA Kit
Reactivity	Human
Range	3.12 pg/mL-200 pg/mL
Target	SFRP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q8N474
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.