

## Product Datasheet

### Human CAPN2 ELISA Kit (orb406277)

**Description** Human calpain 2, (m/II) large subunit (CAPN2) ELISA kit

**Reactivity** Human

**Range** 0.156 ng/mL-10 ng/mL

**Target** CAPN2

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, cell lysates

**Uniprot ID** **P17655**

**Sensitivity** 0.039 ng/mL

**Expiration Date** Please enquire.