

Product Datasheet

Human MASP1 ELISA Kit (orb406293)

Description	Human mannan-binding lectin serine peptidase 1 (C4/C2 activating component of Ra-reactive factor) (M...
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	MASP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P48740
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.