

## Product Datasheet

### Human GAPDH/G3PDH ELISA Kit (orb406296)

|                          |                                                                        |
|--------------------------|------------------------------------------------------------------------|
| <b>Description</b>       | Human glyceraldehyde-3-phosphate dehydrogenase (GAPDH/G3PDH) ELISA kit |
| <b>Reactivity</b>        | Human                                                                  |
| <b>Range</b>             | 31.25 pg/mL-2000 pg/mL                                                 |
| <b>Target</b>            | GAPDH/G3PDH                                                            |
| <b>Note</b>              | For research use only                                                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates                        |
| <b>Uniprot ID</b>        | <b>P04406</b>                                                          |
| <b>Sensitivity</b>       | 7.81 pg/mL                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                        |