

Product Datasheet

Human HBEGF ELISA Kit (orb406305)

Description	Human heparin-binding epidermal growth factor-like growth factor,HB-EGF ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	HBEGF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q99075
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.