

Product Datasheet

Human BIRC2 ELISA Kit (orb406312)

Description	Human Baculoviral IAP repeat-containing protein 2(BIRC2) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Target	BIRC2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13490
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.