

Product Datasheet

Human BIRC3 ELISA Kit (orb406313)

Description	Human Baculoviral IAP repeat-containing protein 3(BIRC3) ELISA kit
Reactivity	Human
Range	37.5 pg/mL-2400 pg/mL
Target	BIRC3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q13489
Sensitivity	9.4 pg/mL
Expiration Date	Please enquire.