

Product Datasheet

Human MX1 ELISA Kit (orb406317)

Description	Human Interferon-induced GTP-binding protein Mx1(MX1) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	MX1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, tissue homogenates, cell lysates
Uniprot ID	P20591
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.