

Product Datasheet

Human HMW Adiponectin ELISA Kit (orb406324)

Description	Human High Molecular Weight Adiponectin(HMW Adiponectin)ELISA Kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	HMW Adiponectin
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, urine
Uniprot ID	Q15848
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.