

Product Datasheet

MFSD6 antibody (orb40633)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MFSD6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88088
Uniprot ID	Q6ZSS7
NCBI	NP_060164.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

