

Product Datasheet

Human GOT2 ELISA Kit (orb406339)

Description	Human glutamic-oxaloacetic transaminase 2, mitochondrial (aspartate aminotransferase 2) (GOT2) ELISA...
Reactivity	Human
Range	0.156 mU/mL-10 mU/mL
Target	GOT2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P00505
Sensitivity	0.039 mU/mL
Expiration Date	Please enquire.