

Product Datasheet

Human BMPR1A ELISA Kit (orb406354)

Description	Human Bone morphogenetic protein receptor 1A,BMPR-1A ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	BMPR1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P36894
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.