

Product Datasheet

Human BMP6 ELISA Kit (orb406360)

Description	Human Bone Morphogenetic Protein 6,BMP-6 ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	BMP6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cell lysates
Uniprot ID	P22004
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.