

Product Datasheet

Human SREBF1 ELISA Kit (orb406363)

Description	Human Sterol regulatory element-binding protein 1(SREBF1) ELISA kit
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	SREBF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P36956
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.