

Product Datasheet

Human TIE1 ELISA Kit (orb406399)

Description	Human tyrosine kinase with immunoglobulin-like and EGF-like domains 1,Tie-1 ELISA kit
Reactivity	Human
Range	0.94 ng/mL-60 ng/mL
Target	TIE1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35590
Sensitivity	0.23 ng/mL
Expiration Date	Please enquire.