

Product Datasheet

Human YBX1 ELISA Kit (orb406409)

Description	Human Nuclease-sensitive element-binding protein 1(YBX1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	YBX1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P67809
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.