

Product Datasheet

Human RNASE2 ELISA Kit (orb406411)

Description	Human ribonuclease, RNase A family, 2 (liver, eosinophil-derived neurotoxin) (RNASE2) ELISA kit
Reactivity	Human
Range	78 pg/mL-5000 pg/mL
Target	RNASE2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P10153
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.