

Product Datasheet

Human NFKB1 ELISA Kit (orb406415)

Description	Human Nuclear factor-kappa B,NF-kappaB ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	NFKB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P19838
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.