

Product Datasheet

Human NFE2L2 ELISA Kit (orb406418)

Description	Human Nuclear factor erythroid 2-related factor 2(NFE2L2) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	NFE2L2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, tissue homogenates, cell lysates
Uniprot ID	Q16236
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.