

Product Datasheet

Human MAGEA3 ELISA Kit (orb406421)

Description	Human melanoma-associated antigen 3 (MAGE-3) ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	MAGEA3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, tissue homogenates
Uniprot ID	P43357
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.