

Product Datasheet

Human ALOX12 ELISA Kit (orb406428)

Description	Human Arachidonate 12-lipoxygenase, 12S-type (ALOX12/LOG12) ELISA kit
Reactivity	Human
Range	1.25 ng/mL-80 ng/mL
Target	ALOX12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P18054
Sensitivity	0.31 ng/mL
Expiration Date	Please enquire.