

## Product Datasheet

### Human XDH ELISA Kit (orb406432)

|                          |                                             |
|--------------------------|---------------------------------------------|
| <b>Description</b>       | Human Xanthine Oxidase(XOD)<br>ELISA Kit    |
| <b>Reactivity</b>        | Human                                       |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL                        |
| <b>Target</b>            | XDH                                         |
| <b>Note</b>              | For research use only                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                |
| <b>Sample Types</b>      | serum, plasma, urine, tissue<br>homogenates |
| <b>Uniprot ID</b>        | <b>P47989</b>                               |
| <b>Sensitivity</b>       | 0.039 ng/mL                                 |
| <b>Expiration Date</b>   | Please enquire.                             |