

Product Datasheet

Human IMPDH2 ELISA Kit (orb406441)

Description	Human Inosine-5'-monophosphate dehydrogenase 2 (IMPDH2) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Target	IMPDH2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P12268
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.