

Product Datasheet

Human MMP8 ELISA Kit (orb406468)

Description	Human Matrix metalloproteinase 8/Neutrophil collagenase,MMP-8 ELISA kit
Reactivity	Human
Range	39 pg/mL-2500 pg/mL
Target	MMP8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P22894
Sensitivity	9.75 pg/mL
Expiration Date	Please enquire.