




Product Datasheet

CD1E antibody (orb40647)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CD1E
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43626
Uniprot ID	P15812
NCBI	NP_001036048.1, NP_001036052.1, NP_001172042.1, NP_001172037.1, NP_001036051.1, NP_001172039.1, NP_001172036.1, NP_001172041.1, NP_001172044.1, NP_001036050.1, NP_001172043.1, NP_112155.2, NP_001036049.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MDA-MB435...



Immunohistochemical staining of paraffin...