

Product Datasheet

Human TIMP3 ELISA Kit (orb406470)

Description	Human tissue inhibitors of metalloproteinase 3,TIMP-3 ELISA Kit
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Target	TIMP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35625
Sensitivity	7.8 pg/ml
Expiration Date	Please enquire.