

Product Datasheet

Human SDF-1 β /SDF1B ELISA Kit (orb406473)

Description	Human stromal cell derived factor 1 β (SDF-1 β /SDF1B) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	SDF-1 β /SDF1B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P48061
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.