

## Product Datasheet

### Human KLK5 ELISA Kit (orb406480)

|                          |                                                                       |
|--------------------------|-----------------------------------------------------------------------|
| <b>Description</b>       | Human Kallikrein-5(KLK5) ELISA kit                                    |
| <b>Reactivity</b>        | Human                                                                 |
| <b>Range</b>             | 62.5 pg/mL-4000 pg/mL                                                 |
| <b>Target</b>            | KLK5                                                                  |
| <b>Note</b>              | For research use only                                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                          |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates,<br>tissue homogenates, saliva |
| <b>Uniprot ID</b>        | <b>Q9Y337</b>                                                         |
| <b>Sensitivity</b>       | 15.6 pg/mL                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                       |