

Product Datasheet

Human CHIT1 ELISA Kit (orb406488)

Description	Human chitinase 1 (chitotriosidase) (CHIT1)ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	CHIT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q13231
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.