

## Product Datasheet

### Human PTH 1-34 ELISA Kit (orb406494)

|                          |                                                    |
|--------------------------|----------------------------------------------------|
| <b>Description</b>       | Human parathyroid hormone 1-34(PTH 1-34) ELISA kit |
| <b>Reactivity</b>        | Human                                              |
| <b>Range</b>             | Request Information                                |
| <b>Target</b>            | PTH 1-34 Ab                                        |
| <b>Note</b>              | For research use only                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                       |
| <b>Sample Types</b>      | serum, plasma                                      |
| <b>Sensitivity</b>       | Request Information                                |
| <b>Expiration Date</b>   | Please enquire.                                    |