

Product Datasheet

Human ITIH1 ELISA Kit (orb406504)

Description	Human inter-alpha-trypsin inhibitor heavy chain H1(ITIH1) ELISA Kit
Reactivity	Human
Range	15.6 ng/mL-1000 ng/mL
Target	ITIH1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, tissue homogenates
Uniprot ID	P19827
Sensitivity	3.9 ng/mL
Expiration Date	Please enquire.