

Product Datasheet

Human CTHRC1 ELISA Kit (orb406511)

Description	Human Collagen triple helix repeat-containing protein 1(CTHRC1) ELISA kit
Reactivity	Human
Range	6.25 ng/mL-400 ng/mL
Target	CTHRC1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q96CG8
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.