

Product Datasheet

Human ADAM12 ELISA Kit (orb406527)

Description	Human A Disintegrin And Metalloprotease 12,ADAM12 ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	ADAM12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	O43184
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.