

Product Datasheet

Human ADAM33 ELISA Kit (orb406530)

Description	Human A Disintegrin And Metalloprotease 33,ADAM33 ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	ADAM33
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, cell culture supernates, saliva, tissue homogenates
Uniprot ID	Q9BZ11
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.