

Product Datasheet

Human ADAM9 ELISA Kit (orb406532)

Description	Human A Disintegrin And Metalloprotease 9,ADAM9 ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	ADAM9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates, ascitic fluid
Uniprot ID	Q13443
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.