

Product Datasheet

Human ADAMTS1 ELISA Kit (orb406534)

Description	Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 1(ADAMTS1) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20ng/mL
Target	ADAMTS1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, urine
Uniprot ID	Q9UHI8
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.