

Product Datasheet

Human ARG2 ELISA Kit (orb406538)

Description	Human Arginase 2(ARG2) ELISA Kit
Reactivity	Human
Range	0.47 mU/mL-30 mU/mL
Target	ARG2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, cerebrospinal fluid (CSF)
Uniprot ID	P78540
Sensitivity	0.12 mU/mL
Expiration Date	Please enquire.