

## Product Datasheet

### Human SSA/Ro ELISA Kit (orb406565)

|                          |                              |
|--------------------------|------------------------------|
| <b>Description</b>       | Human SSA/Ro ELISA Kit       |
| <b>Reactivity</b>        | Human                        |
| <b>Range</b>             | Request Information          |
| <b>Target</b>            | SSA/Ro Ab                    |
| <b>Note</b>              | For research use only        |
| <b>Application notes</b> | Detection Wavelength: 450 nm |
| <b>Sample Types</b>      | serum, plasma                |
| <b>Sensitivity</b>       | Request Information          |
| <b>Expiration Date</b>   | Please enquire.              |