

Product Datasheet

Human AVP ELISA Kit (orb406612)

Description	Human antidiuretic hormone/vasopressin/arginine vasopressin,ADH/VP/AVP ELISA Kit
Reactivity	Human
Range	1.25 pg/mL-80 pg/mL
Target	AVP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P01185
Sensitivity	0.312 pg/mL
Expiration Date	Please enquire.