

Product Datasheet

Human AMHR2 ELISA Kit (orb406634)

Description	Human Anti-Muellerian hormone type-2 receptor(AMHR2) ELISA kit
Reactivity	Human
Range	0.32 ng/mL-12 ng/mL
Target	AMHR2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q16671
Sensitivity	8 pg/mL
Expiration Date	Please enquire.