

## Product Datasheet

### Human smooth muscle ELISA Kit (orb406643)

|                          |                                    |
|--------------------------|------------------------------------|
| <b>Description</b>       | Human smooth muscle,ASMA ELISA Kit |
| <b>Reactivity</b>        | Human                              |
| <b>Range</b>             | Request Information                |
| <b>Target</b>            | ASMA                               |
| <b>Note</b>              | For research use only              |
| <b>Application notes</b> | Detection Wavelength: 450 nm       |
| <b>Sensitivity</b>       | Request Information                |
| <b>Expiration Date</b>   | Please enquire.                    |